

Rapid Design and Manufacture of Single Species, Site-Specific ADCs using Cell-Free Synthesis with nnAAs

From Concept to Clinic in 365 days ?

Trevor J. Hallam, PhD Chief Scientific Officer 3rd World ADC Summit, San Francisco October 24th, 2012

ADCs – Exciting prospects but . . .





Heterogeneity translates to poor PK, stability and efficacy



The Non-Natural Amino Acid Advantage

- Controlled stability: nnAA chemical space provides alternatives to cysteine or lysine for creating stable MAb~drug junction
- 2. Homogeneity: site-specific conjugation using orthogonal chemistries regulates number and location of drugs attached to Mab



Anatomy of a Disruptive Protein Design and Manufacturing Solution....

- Rapid Make-Test
 - Fast protein production (from DNA to g's/L in hours)
 - a 1-week cycle time for purification, characterization and testing
 - Micro-titre plate evaluation of 100's of variants in parallel
 - Parallel evaluation of optimal production conditions
- Chemical Diversity
 - Incorporation of non-natural amino acids
- Rapid Scale-up
 - Linearly scalable and predictable expression allows for seamless and rapid progression to specificity testing, efficacy models and GLP tox studies
- Same system for cGMP manufacture







 Examples of Ab Fragments Produced by CF synthesis

Typical Expression Levels:

- 0.25 -1 Gram / Liter (Pre-Optimized)
- 6-8 Hour Reaction
- 30% Extract

Translation of nnAA-Containing Proteins : Enables Site-Specific Conjugation

Su

l **t r O** Biopharma

Data driven design: production of many variants in hours

Surface Scan

- Mutate Sites in IgG: Choose nnAA sites using rational design, or just make all of them!
- Produce nnAA IgG: Incorporate nnAA at 100's of chosen sites
- Conjugate: Conjugate nnAA with appropriate chemistry
- Purify: Separate conjugated IgG away from unincorporated linkerwarhead
- Test: Assay conjugated IgG's for binding and cell killing

TAG Screen

Cu Free Click Conjugation Chemistry

Rapid Selection of Optimal Sites for Expression, Conjugation, Binding and Killing

Sutro

Site-Specific Incorporation of MMAF Yields Potent Cytotoxic ADC

- Comparable IC50 to published ADC at equivalent IgG concentration
- More potent at equivalent drug concentration

Pharmacokinetics of Cell Free Produced ADCs are Comparable to Cell-Based derived ADCs

CF-Trastuzumab Drug conjugate pharmacokinetics are in good agreement with Trastuzumab- MMAF conjugate literature

	Sutro Data ^a	SeaGen Data ^b
	CF-Trastuzumab Drug Conjugate	Trastuzumab MMAF Conjugates
AUC _{inf} [day/µg/mL]	248	299
Clearance [mL/d/kg]	8.1	9
Half-life [d]	8	10

^a5 mg/kg, Balb-c mice, Sutro ^b2mg/kg, Sprague-Dawley rats (US7994135B2, SeaGen patent 2011)`

Efficacy Study: Agly Trastuzumab DC Dose Response

Site Dependent Impact on Cell Killing Observed

TAG Scanning: DAR vs. Cell Killing

Efficient Integration of nnAAs and Conjugation Chemistry are key for manufacturing ADCs

- nnAA Incorporation Efficiency
 - Site Dependence
 - RF1 Attenuation
 - Multiple nnAA Incorporation
- Conjugation Efficiency
 - Site Dependence
 - Improving Conjugation Kinetics

nnAA Incorporation Efficiency; Site Dependence in HC

nnAA Incorporation Efficiency; Site Dependence in LC

nnAA Incorporation Efficiency:

RF1 Attenuation

- nnAA Incorporation Efficiency
- The second second

Protease Sensitive RF1 is inactivated during extract production

Abbrev.	Sample	
М	Marker	
Р	Pellet	
L	Lysate	
С	Clarified	
1	1 hour time point	
2	2 hour time point	
3	3 hour time point	

De-compartmentalized Extract Cleaves RF1

New Strains Have Similar Growth Rates To: Control Strain

High Growth Rate is Historically Predictive of Good Extract Performance

RF1⁻ Strain Engineered Extract Boosts nnAA-Protein Production

Strain: RF1 Attenuating Modification

Suppression of amber codon in HC

Incorporation of p-Acetyl Phe into IgG at three selected sites

Su

l **t r O** Biopharma

RF1 Depletion through Strain Engineering: Sutro Extract Can Produce IgG w/ Multiple nnAA's Sutro

- Incorporation of multiple nnAAs in each IgG LC and HC
- Enables Multiple (4,6,8,10+) <u>site-specific</u> combinations of warheads/lgG
- Intractable using cell-based systems or wildtype extract due to significant accrued losses in yield
- Combinations of sites screened for:
 - nnAA Incorporation efficiency/expression
 - Conjugation Efficiency
 - Stability
 - PK/Potency in vivo

4 nnAA IgG Yields Using Combinations of 2 LC and 2 HC Sites

- nnAA: nn-AA1
- Synthetase 3

- nnAA: pN3nnAA
- Synthetase 1

Several 4 nnAA combos (2LC+2HC) express similarly to WT

Conjugation Efficiency

Cu Free Click Conjugation Chemistry

ō, ö

ö

0

ČO₂H

Conjugation Efficiency

Synthesis of site-specific antibody-drug conjugates using unnatural amino acids

Jun Y. Axup^{a,b}, Krishna M. Bajjuri^c, Melissa Ritland^d, Benjamin M. Hutchins^{a,b}, Chan Hyuk Kim^{a,b}, Stephanie A. Kazane^{a,b}, Rajkumar Halder^{a,b}, Jane S. Forsyth^d, Antonio F. Santidrian^d, Karin Stafin^d, Yingchun Lu^e, Hon Tran^e, Aaron J. Seller^e, Sandra L. Biroc^e, Aga Szydlik^e, Jason K. Pinkstaff^e, Feng Tian^e, Subhash C. Sinha^c, Brunhilde Felding-Habermann^d, Vaughn V. Smider^{c,1}, and Peter G. Schultz^{a,b,1}

^aDepartment of Chemistry, ^cDepartment of Molecular Biology, ^dDepartment of Molecular and Experimental Medicine, and ^bThe Skaggs Institute for Chemical Biology, The Scripps Research Institute, La Jolla, CA 92037; and ^eAmbrx, Inc., La Jolla, CA 92037

Contributed by Peter G. Schultz, July 20, 2012 (sent for review January 27, 2012)

Conjugation. Oxime ligation between pAcPhe on the antibody and alkoxyamine-functionalized auristatin was carried out in 100 mM acetate buffer pH 4.5 with 100 μM (5 mg/mL) Fab and 3 mM AF (30-fold excess) for 1–2 d at 37 °C. **IgG was conjugated at 66.7 μM (10 mg/mL) and 1.3 mM AF** (20-fold excess) for 4 d at 37 °C.

Novel nnAAs with Boosted Conjugation Kinetics

New Chemistry offers Improved Kinetics, Flexibility

Some sites are completely conjugated in under 4 hours!

Cell-Free Manufacturing: GMP Facility a Critical Element

Rapid Production of Biotherapeutics

Su

tro

- Direct linear scale-up from HTS to production scale
- Uses standard bioreactors & downstream equipment
- Minimal, rapid process development
- Gene sequence to drug substance in days

Sutro's technology enables Best-in-Class ADCs....

- Data-driven assessment of many potential ADC variants ensures optimized positioning for suppression, expression, conjugation efficiency, payload, cell killing, stability and pharmacokinetic properties
- Flexible and rapid scaling means pharmacodynamic and exploratory toxicology assessment studies can be front-loaded into the discovery phase
- Material for GLP Tox studies and clinical studies can be delivered in days/weeks from selection of Clinical Development Candidate
- 365 days?